



**FlowBoot**<sup>™</sup>  
Memory Soft Compression Technology

# ***Cold Compression FlowBoot***<sup>™\*</sup>

\*U.S. and Foreign Patent Pending  
**MADE IN USA**

- ✓ New Memory Soft Compression Technology
- ✓ No constricting hoses or cables
- ✓ No annoying readjusting
- ✓ Long lasting Cold Compression
- ✓ Easy to use design

956-551-2828 • [FlowBoot.com](http://FlowBoot.com)

# Cold Compression FlowBoot™\*


## Directions

1. Soak boots in water over night, then place individually in a plastic bag.
2. Stuff the top of the bag inside the middle portion of the boots to keep sides from touching and place in freezer until frozen.  
*(The plastic keeps the boots from sticking to each other and to the freezer.)*
3. When ready to use, thaw the boot slightly in water or add a little water to plastic.
4. Once pliable, strap on boots with the wrap.
5. Leave the boots on 20 to 30 minutes.  
*(20 minutes is what the medical field recommends, but if there is an injury or hot spot, reapply multiple times.)*
6. If on the road, keep the boots in ice water.

\*It's best to use when the horse's body is already cooled off. Rubbing alcohol can be used to disinfect and will help the cold. Water is the conducting element and transfers the cold, best when wet.



- Cold compression is significantly colder than ice alone due to improved skin contact.
- Tissue reaches its lowest temperature faster and the tissue maintains its coolness even after treatment ends.
- Flowboot targets the superficial digital flexor tendon and proximal sesamoid area where most problems occur.
- Keeping the heat and swelling out promotes recovery and prevents soreness.



**Our Guarantee:** We take great pride in the quality of our products and stand behind everything we sell. If you are not completely satisfied with your purchase, you can return it for full product refund. Items must be returned within 45 days of purchase.  
**Warning:** Wraps are not to be used without the boots. This can cause damage to tendons.

**FOR QUESTIONS, CONTACT: [WWW.FLOWBOOT.COM](http://WWW.FLOWBOOT.COM) OR CALL 956.551.2828**



# Cold Compression FlowBoot<sup>TM\*</sup>


## Directions

1. Soak boots in water over night, then place individually in a plastic bag.
2. Stuff the top of the bag inside the middle portion of the boots to keep sides from touching and place in freezer until frozen.  
*(The plastic keeps the boots from sticking to each other and to the freezer.)*
3. When ready to use, thaw the boot slightly in water or add a little water to plastic.
4. Once pliable, strap on boots with the wrap.
5. Leave the boots on 20 to 30 minutes. *(20 minutes is what the medical field recommends, but if there is an injury or hot spot, reapply multiple times.)*
6. If on the road, keep the boots in ice water.

\*It's best to use when the horses body is already cooled off.

\*Rubbing alcohol can be used to disinfect and will help the cold.

\*Water is the conducting element and transfers the cold, best when wet.



WARNING: WRAPS ARE NOT TO BE USED WITHOUT THE BOOTS. THIS CAN CAUSE DAMAGE TO TENDONS.

FOR QUESTIONS, CONTACT: [WWW.FLOWBOOT.COM](http://WWW.FLOWBOOT.COM) OR CALL 956.551.2828